CitoVac

Instruction Manual

Manual No.: 15927001

Date of Release: 01-01-2021

Instruction Manual

Table of Contents	Page
User's Guide	1
Reference Guide.....	23

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the Serial No. and Voltage on the type plate of the machine itself. We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

Instruction Manuals: Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals: Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to changes without notice. The manual may mention accessories or parts not included in the present version of the equipment.

Original instructions. The contents of this manual is the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2011 .

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup
Denmark
Telephone +45 44 600 800
Fax +45 44 600 801

CitoVac

Safety Precaution Sheet


To be read carefully before use

1. The operator should be fully aware of the use of the units according to the Instruction Manual. The unit must be placed in an adequate working position.
2. Never use vacuum lids other than the ones delivered from Struers.
3. If any cracks occur in the vacuum lid it must be replaced immediately.
4. Ensure that the compressed air and vacuum are switched off before removing the hoses.
5. Use only consumables which are suitable for vacuum impregnation. Mounting materials must be handled with care. See the Instruction Manual and Material Safety Data Sheet for the resin.

The equipment should only be used for its intended purpose and as detailed in the Instruction Manual.

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.).

Disposal

Equipment marked with a WEEE symbol  contain electrical and electronic components and must not be disposed of as general waste.

Please contact your local authorities for information on the correct method of disposal in accordance with national legislation.

User's Guide

Table of Contents	Page
1. Getting Started	
Checking the Contents of Packing	2
Consumables	2
Placing CitoVac.....	2
Getting Acquainted with CitoVac.....	3
Supplying Power	4
Connection to the Machine.....	4
Connection to Compressed Air/ Vacuum.....	5
CitoVac for External Vacuum Pump.....	5
2. Operation	
Using the Controls.....	7
Control Panel of CitoVac.....	7
Control Panel Functions	7
Main switch	7
Software Settings.....	8
Using Software Menus	9
Configuration Menu.....	10
Process Configuration.....	10
Options	11
Setting the Language.....	12
Selecting a Method	13
Editing a Method	14
Preparing the Impregnation.....	15
Checking the Vacuum Chamber	16
Impregnation	18
Running a Process.....	20
Gluing	21

1. Getting Started

Checking the Contents of Packing

In the packing box you should find the following parts:


- 1 CitoVac unit
- 2 Power cables
- 1 Universal mounting cup holder (inside vacuum chamber)
- 1 Hose for vacuum
- 1 Coupling for vacuum hose
- 2 Hose clamps
- 1 Washer (sealing ring for vacuum chamber)
- 1 Instruction Manual set

Consumables

- 1 Chamber protector (inside vacuum chamber)
- 1 Mixing cup
- 1 Dispensing tube
- 1 Clip for dispensing tube

Placing CitoVac

- CitoVac should be placed on a sturdy table or work bench at an appropriate working height.
- Level CitoVac by means of the adjustable leg.

- Place CitoVac close to connections for mains power.
- Place CitoVac close to the outlet for compressed air/ vacuum.

NB

Impregnation materials such as epoxy may produce fumes. CitoVac should only be operated where ample ventilation is possible, e.g. in a fume cupboard.

**Getting Acquainted
with CitoVac**

Take a moment to familiarise yourself with the location and names of the CitoVac components.

- ① Pivot joint of the lid
- ② Lid
- ③ Washer
- ④ Vacuum chamber
- ⑤ Wheel for rotating mounting cups
- ⑥ Front Panel
- ⑦ Vacuum valve
- ⑧ Dispensing tube
- ⑨ Mixing cup holder

Supplying Power

Always remember to switch the power off when installing electrical equipment.

DANGER!


The machine must be earthed

IMPORTANT

Check that the mains voltage corresponds to the voltage stated on the type plate on the back of the machine, 100-240V 50-60Hz.

CitoVac is shipped with 2 types of Mains cables:

European plug


The 2-pin (European Schuko) plug is for use on single-phase connections.

If the plug supplied on this cable is not approved in your country, then the plug must be replaced with an approved plug. The leads must be connected as follows:

Yellow/green: earth
Brown: phase (live)
Blue: neutral

American plug

The 3-pin (North American NEMA) plug is for use on single-phase power connections.

If the plug supplied on this cable is not approved in your country, then the plug must be replaced with an approved plug. The leads must be connected as follows:

Green: earth
Black: phase (live)
White: neutral

Connection to the Machine

Both cables are on the other end equipped with an IEC 320 cable connector that has to be connected to CitoVac.

Connection to Compressed Air/ Vacuum

- Mount the coupling on the vacuum hose and secure with a hose nipple.
- Connect the vacuum hose to the rear of CitoVac.
- Mount the hose on the compressed air/ vacuum supply.

CitoVac for External Vacuum Pump

CitoVac models for use with an external vacuum pump can be connected to the vacuum pump using a Cooli-1 Control Unit.

- Connect the CitoVac and the Cooli-1 Control Unit using the communication lead supplied.

- ① Fuse
- ② Main switch
- ③ Power supply
- ④ Communication cable to Control Unit
- ⑤ Communication socket
- ⑥ Compressed air/ Vacuum connection

Rear of Cooli Control Unit

- Mount a four pole plug on the power lead from the vacuum pump and connect to the Cooli-1 Control Unit.
NB Four pole plug must be ordered separately, Cat. No: 2XM10221. The plug must be assembled by a qualified technician.

- ⑦ Power supply to Vacuum pump (four pole plug)
- ⑧ Mains power cable
- ④ Communication cable, connected to CitoVac

Vacuum / Compressed Air Hose

Take care that the hose does not have any sharp bends. If a longer hose is required it is recommended to choose a hose with a larger internal diameter, in order to secure a sufficient flow.

2. Operation

Using the Controls Control Panel of CitoVac

Control Panel Functions

Name	Key	Function	Name	Key	Function
UP MENU KEY		Moves highlighted menu item upwards and increases values of selected parameters, and moves cursor left in menus.	DOWN MENU KEY		Moves highlighted menu item downwards and decreases values of selected parameters and moves cursor right in menus.
ESCAPE	Esc	Moves one step backward in menus.	ENTER		Selects highlighted menu items and enters (saves) modified parameter values.
START		Starts the vacuum process.	STOP		Stops the vacuum process: Press once to pause the process, press twice to stop the process.

Main switch

The Main switch is located at the rear of CitoVac.

Software Settings

Initial Start Up Screen


When CitoVac is initially powered up, using the mains switch, the following two screens will be displayed:

Please Note

The sample screens in this Instruction Manual show a number of possible texts. The actual display screen may differ from the samples in the Manual.


A pop-up will appear to select the preferred Language setting:

Use the up  and down  menu keys to highlight the required language then press the  enter key to select it.

Main Menu

The MAIN MENU will automatically replace these screens.

Using Software Menus

Use the up  and down  menu keys, the escape Esc and enter  keys to navigate and use CitoVac's software menus

Up and Down Menu Keys

- Press the up and down menu keys to highlight menu items.

ESC Key

Esc

- Press the ESC key to return to the previous menu.
- Press the ESC key to exit a selected item.

Enter Key

- Press the enter key to select a highlighted item.
- Press the enter key to save a value that has been changed

Configuration Menu

Two sub menus that can be accessed from the CONFIGURATION menu *Process Configuration* and *Options*.

Process Configuration

PROCESS CONFIG	
Units	: Bar
Ambient pressure	: 1.02
Pressure mode	: Relative
Keep vacuum	: Off

Units

Units can be set to Bar (default), kPa or PSI.

Ambient Pressure

CitoVac measures the difference in relative pressure and not the actual absolute pressure.

The actual value of ambient pressure can be input for greater accuracy. This is particularly important when using *Absolute Pressure* mode (see below).

Pressure Mode

Pressure mode can be set to *Absolute* or *Relative*.

The difference in pressure between ambient pressure and the pressure inside the chamber is measured.

Relative: *Vacuum* is shown in the method display.

METHODS		
	Vacuum (Bar)	Time (h:m:s)
Method A	0.60	1:00:00
Method B	0.90	0:30:00
Method C	0.40	0:10:00

Absolute: Absolute pressure is calculated using the value of *Ambient pressure* set in PROCESS CONFIG. *Pressure* is shown in the method display.

METHODS		
	Pressure (Bar)	Time (h:m:s)
Method A	0.42	1:00:00
Method B	0.12	0:30:00
Method C	0.62	0:10:00

Keep Vacuum

CitoVac can be set to *On* to maintain vacuum after the process is complete or to *Off* to release the vacuum. The default is *Off*.

Options

OPTIONS	
Language	: English
Display contrast	: 60
Acoustic signal	: Off

Language

The language can be set to English (default), German, French, Spanish or Italian.

Display contrast

The contrast settings of the display can be adjusted to suit individual preferences (default value: 60, adjustment range: 1-100).

Acoustic signal


The option to have the machine signal audibly when it has completed the mounting process for a sample can be switched **On** or **Off**. If switched **On**, CitoVac will “beep” when keypads are pressed.

Setting the Language

When the machine is first started the language selected will be English.

To change the language after the initial start up:


From the *MAIN MENU*, use down menu key to highlight *Config*.

Push the enter key to activate the *CONFIGURATION* menu.


Use the down menu key to highlight *Options* then press the enter key to activate the *OPTIONS* menu.

Push the enter key to activate the Language pop-up.

Use the up  and  down menu keys to highlight the required language then press the  enter key to select it.


Press Esc key to return to the *MAIN MENU*.

Selecting a Method

CitoVac has three built-in methods where the pressure/ vacuum and process time can easily be set and saved.

From the *MAIN MENU*, use  down menu key to highlight *Methods*.

Push the  enter key to activate the *METHODS* menu. The *Vacuum* and *Process Time* settings for each method will be shown.

METHODS		
	Vacuum (Bar)	Time (h:m:s)
Method A	0.60	1:00:00
Method B	0.90	0:30:00
Method C	0.40	0:10:00

Editing a Method

Each of the three methods can be edited and saved.
From the *METHODS* menu:

Use the up and down menu keys to highlight the required Method.

Push the ↵ enter key to open the Method window.

Method A – Ready		
	Pressure (Bar)	Time (hh:mm:ss)
Set	0.42	1:00:00
Actual	✓	1:00:00

Use the up and down menu keys to highlight the setting to be edited e.g. *Vacuum*.

Push the ↵ enter key to open the editing scroll bars.

↓

Method A – Ready		
	Pressure (Bar)	Time (hh:mm:ss)
Set	0.42	1:00:00
Actual	✓	1:00:00

Use the up and down menu keys to adjust the value.

Push the ↵ enter key to accept the new value.

Preparing the Impregnation

- Place the clean and dry specimen in a suitable mounting cup.
NB Before impregnation the specimens should be cleaned and degreased.
- Check that the chamber protector is in place.
- Place the mounting cup in the vacuum chamber.
- Check that each of the cups will be directly under the nozzle of the dispensing tube by turning the rotating wheel.


Note:

Struers have a holder specially designed for Struers mounting cups. The mounting cup holder is delivered with a support for safe storage when not placed in the vacuum chamber.

Checking the Vacuum Chamber

- Before operation, check that the lid is not cracked or has fissures, or it might implode when exposed to vacuum.
- Check that the washer is clean and undamaged.
- Put the nozzle of the dispensing tube through the opening in the vacuum chamber and press firmly into place.

- Position the dispensing tube in the groove of the Vacuum valve.
 - **Tip:** The valve must be fully open (the line on the valve should be facing front). Stretch the tube slightly to ease correct positioning in the groove.

- Turn the lid so that it is directly over the chamber.
- Close the valve and press Start.

- Press down on the pivot joint of the lid until it forms a seal with the vacuum chamber.


Ensure that the lid is directly over the chamber and that there are no leaks around the edges. Should some leaking occur, release the vacuum, reposition the lid and reapply vacuum.

Impregnation

For mounting, Struers manufactures a number of mounting cups that can be used without releasing agent.

- When the specimens have been under vacuum for an appropriate time - from a couple of minutes for not very porous specimens up to half an hour for very porous specimens - the impregnation can begin:
- Place the Mixing Cup with the ready-mixed impregnation material in the holder.
- Place the end of the dispensing tube in the bottom of the cup and secure with the clip placed on the lip of the mixing cup. Check that the tube is near the bottom of the cup, or air may be sucked into the system and create splashing in the vacuum chamber.

- Turn the Holder with the mounting cups until one of the cups is directly under the dispensing tap.
- Slowly open the vacuum valve until the impregnation material runs into the cup with an appropriate flow.

- When the specimen is covered with impregnation material, turn off the flow, turn the holder so that the next mounting cup is under the dispensing tap and repeat the filling procedure. Repeat until all mounting cups have been filled up. If necessary, refill the mixing cup with impregnation material.

Tip:

Some impregnation material, e.g. Struers EpoFix may be heated a little (max. 40°C/104°F) before pouring it over the specimen. This gives the mixture a lower viscosity, thus ensuring a more thorough penetration into the pores of the specimen.

Running a Process

From the METHODS menu:

Use the up and down menu keys to highlight the required Method.

Push the enter key to open the Method window.

Method A – Ready		
	Pressure (Bar)	Time (hh:mm:ss)
Set	0.42	1:00:00
Actual	✓	1:00:00

Push Start to start the process.

Method A – Running		
	Pressure (Bar)	Time (hh:mm:ss)
Set	0.42	1:00:00
Actual	✓	0:58:30
100%		

Press stop once, to pause the process.

Method A – Paused		
	Pressure (Bar)	Time (hh:mm:ss)
Set	0.42	1:00:00
Actual	✓	0:58:29
Press start to resume		


Method settings can be edited whilst the process is paused.

E.g. if the vacuum is too high and the resin is seen to start boiling whilst developing a method, the vacuum setting can be reduced whilst running the process.

- ⊖ Press Stop again to stop the process.

Method A – Stopped		
	Pressure (Bar)	Time (hh:mm:ss)
Set	0.42	1:00:00
Actual	--	1:00:00
New proc.: Push Start		

Gluing Accessory

- ① Lid
- ② Support ring
- ③ Pressing rod
- ④ Spare rubber plug for pressing rod

- Place the lid with the pressing rod on CitoVac.
- Position the support ring underneath the holder so that it does not move when gluing.
- Apply a suitable amount of resin to the specimen. One drop is usually enough.
- Place the specimen with the glued surface up, and place a glass slide on the specimen.
- Set vacuum to maximum (pressure to minimum) and start CitoVac.
- Press the glass slide with the pressing rod and gently move the slide back and forth. When the specimen is firmly stuck to the slide, stop CitoVac.
- Remove the slide with glued specimen and leave to harden.

Reference Guide

Table of Contents	Page
1. Advanced Operation	
Hints	23
Maintenance.....	23
Cleaning the Lid	23
Changing the Washer	23
Service Menu	24
Trouble Shooting.....	24
Vacuum Quality Test.....	24
Consumables and Accessories	25
Technical Data	26
Air supply	26
Vacuum Connection.....	26
Dimensions	26

1. Advanced Operation

Hints

At too high vacuum some of the components in the epoxy evaporate and it looks as if it boils. This can be avoided by reducing the vacuum.

The impregnation process can be paused whilst the vacuum setting is adjusted accordingly.

Maintenance

Cleaning the Lid

Clean the lid periodically with ethyl alcohol. Do NOT use Acetone as this may result in cracks in the lid.

Changing the Washer


Check the vacuum chamber washer at regular intervals to ensure there is no wear or damage. Should there be problems with persistent leaks, this may indicate that the washer needs replacing.

To insert a new washer:

- Hold the two ends of the new washer together and slot into the groove in the vacuum chamber.
- Carefully press the washer into the groove making sure that it is seated evenly around the chamber.
- Check the vacuum to ensure there are no leaks. Should leaks occur, take the washer out and reposition.

Service Menu

The Service menu is accessed from the Main Menu.

Vacuum quality test: To test the vacuum

Reset configuration: To reset configuration to factory settings

Reset Methods: To reset Methods to factory settings

Trouble Shooting

Whilst running a process the screen will display whether or not the set vacuum has been reached.


✓: Vacuum setting has been reached

*: Vacuum setting has not been reached

A message will appear to inform the User that the vacuum setting has not been reached and will allow the User to choose to proceed with the process or to stop.

Vacuum Quality Test

To check the vacuum, perform a Vacuum quality test.

- Check that the vacuum valve is closed properly.
- Check that the nozzle of the dispensing tube is firmly in place.
- Check the washer in the vacuum chamber.
- Check that the compressed air/ vacuum supply is functioning adequately, for example, that there are no bends or kinks in the hose.

If the problem persists, contact a Struers Service Technician.

Consumables and Accessories Please refer to the [Struers Cold Mounting Accessories brochure](#) and the [CitoVac brochure](#) for details of the range available.

Specification	Cat. No.
<i>Cooli-1 Control Unit</i> For controlling of external vacuum pump four pole plug to connect the lead from the pump to Cooli 1.	05761116 2XM10221

Technical Data

Air supply

Compressed air: 4.5-6 bar
Air quality: The air supplied must be of Class-3.6.4 or better, as specified in ISO 8573-1
Consumption of compressed air: 12.5 l/min
Hose connection: $\varnothing 1/4$ "
Vacuum (at compressed air 6 bar): -860 mbar (-655 mm Hg)

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
Date of Release 05.03.2014

CitoVac
Spare Parts and Diagrams

**Always state *Serial No* and *Voltage/frequency*
if you have technical questions or when ordering spare parts.**

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Instruction Manuals: Struers Instruction Manuals may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals: Struers Service Manuals may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to change without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual are the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2014.

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup
Denmark
Telephone +45 44 600 800
Fax +45 44 600 801

Spare Parts and Diagrams

Table of contents	Drawing
<i>Drawings</i>	
CitoVac Complete	15920001F
Cabinet complete	15920010E
Base	15920011D
Front plate complete	15920012C
Lid - complete.....	15920013B
Valve Assembly.....	15920014B
 <i>Diagrams</i>	
Air diagram, CitoVac with ejector	15922000B
Air diagram, CitoVac without ejector	15922001B
Block Diagram, CitoVac with ejector	15923050A
Circuit Diagram, CitoVac with ejector.....	15923100A
Wiring Diagram, CitoVac with ejector.....	15923450B


Some of the drawings may contain position numbers
not used in connection with this manual.

Spare Part list for CitoVac

Drawing	Pos.	Spare Part	Cat no.
15920001		CitoVac Complete	
	10	Cabinet complete	15920010
	20	Base	15920011
	30	Front plate complete	15920012
	40	Lid - complete	15920013
	60	Rubber foot tube	15920102
	80	Cup holder - painted	15929194
	240	Valve - assembly	15920014
	300	Sealing profile - 5x10 - white	2IP25009
	330	Coupler dir 4-M5, KQ2H04-M5	2NF10009
	340	Coupler dir 6-M5, KQ2H06-M5	2NF10010
	350	PU hose $\varnothing 2.5/\varnothing 4$ TU 0425 BU (blue)	2NU14425
15920010		Cabinet complete	
	10	Cabinet - painted	15929150
	70	MOTOR SIDE DAMPER	2GB00010
15920011		Base	
	30	Mains socket with fuse, switch	2XN32107
	50	Vac. ejector EZM13-2-S-A-H-F-Q	2YP10002
	60	SOLENOID VALVE. VDW31,24V DC 1/4	2YM10034
	80	MOTOR SIDE DAMPER	2GB00010
	260	PU Hose $\varnothing 5/\varnothing 8$ TU 0805BU (blue)	2NU14805
	270	PU Hose $\varnothing 5/\varnothing 8$ TU 0805BU (blue)	2NU14805
	280	PU Hose $\varnothing 5/\varnothing 8$ TU 0805BU (blue)	2NU14805
	300	PU Hose $\varnothing 5/\varnothing 8$ TU 0805BU (blue)	2NU14805
	330	2.00AT Glass Fuse 6.3x32 250V CSA	2FU13950
15920012		Front plate complete	
	20	Keyboard Foil, CitoVac	R5920174
	70	PCB for CitoVac, tested	15923000

Spare Part list for CitoVac


15920013		Lid - complete	
	30	Lid - handle	15920131
15920015		Lid for gluing	
		Lid complete	15920015
		Cork with ø4 hole	2GT50008
		Pressing rod	14350080
		Support ring	15920021
		Other Spare Parts	
		Universal mounting cup holder (Holder-unit)	R5920513



F	2012.09.14	Loctite 577 added	SPE	2012.09.14	
A	25.04.2007		PPL	25.04.2007	JFR
Revision	Crea. date dd-mm-yy	Revision description	Draw. Init	Appr. date dd-mm-yy	Appr. Init
		Material:	Scale: 1:2.5	Format: A3	Tolerance: DS/ISO 2768- mK
		ID:	Description:		Surface treat.: None
		15920001 Citovac - complet			Rev: F

Pederstrupvej 84
DK-2750 Ballerup/Copenhagen
Denmark
Phone: +45 44 600 800
Fax: +45 44 600 804

2 (can be adjusted ± 2 mm)



SECTION B-B

110
A little High vacuum grease are to be added.





DETAIL E
SCALE 1 : 1

SECTION A-A

E	05-09-08	Pos 110 added.	MPE		
A	27-03-08		PPL	01-11-06	PPL
Revision	Crea. date dd-mm-yy	Revision description	Draw. Init	Appr. date dd-mm-yy	Appr. Init
		Material:	Scale: 1:2.5	Format: A3	Tolerance: DS/ISO 2768- mK Surface treat.: None
ID:	Description: 15920010 Cabinet-complet				Rev: E




Struers
Pedestergvej 84
DK-2750 Ballerup/Copenhagen
Denmark
Phone: +45 44 600 800
Fax: +45 44 600 804

Minimal bend radius of $\varnothing 8$ pneumatic hoses are R20mm


D	2013-12-04	Y1 & Y2 changed to new version	OCR	4.12.2013	JTV
A	25.04.2007		PPL	25.04.2007	JF
Revision	Crea. date yyyy-mm-dd	Revision description	Draw. Init	Appr. date yyyy-mm-dd	Appr. Init
		Material:	Scale: 1:2	Format: A3	Tolerance: DS/ISO 2768- mK
		ID:	Description:		Weight : g
		15920011 base			Rev: D

Pederstrupvej 84
DK-2750 Ballerup/Copenhagen
Denmark
Phone: +45 44 600 800
Fax: +45 44 600 804

SECTION D-D

See wiring set 15923590
Se leadningssaet 15923590

VIEW C


100 15923450 Wiring diagram

Logitech version:

(70) Install special software with Logitech logo.




C	2011-01-04	Logitech version added.	BRY	2011-01-04	JTV
A	2006-11-20		PPL	2006-11-20	PPL
Revision	Crea. date yyyy-mm-dd	Revision description	Draw. Init	Appr. date yyyy-mm-dd	Appr. Init
		Material:	Scale: 1:2	Format: A3	Tolerance: DS/ISO 2768- mK
ID:	Description:				Weight : g
15920012 frontplate-complet					Rev: C

Pederstrupvej 84
DK-2750 Ballerup/Copenhagen
Denmark
Phone: +45 44 600 800
Fax: +45 44 600 804


C					
B	06.11.2007	Add geometric tolerance (E1+D3-4)	PP	06.11.2007	JF
A	15.11.2006		PP	15.11.2006	PP
Rev	Crea. date dd-mm-yy	Revision description	Draw. Init	Appr. date dd-mm-yy	Appr. Init
F	Struers Pederstrupvej 84 DK-2750 Ballerup Copenhagen Denmark Phone: +45 44600 800 Fax: +45 44600 804		Material: Scale: 1:2	Format: A4	Tolerance: DS/ISO 2768 -mK Surface treat.:
	ID: Description: 15920013 Lid-complet				Rev: B

Pos.20 Skrues i bund, og løsnes 1/2 omgang
 Pos.20 to be fastened and turned ccw 180°

20 Omnifit 50M

C					
B	06.11.2007	Omnifit 50M added	SPE	06.11.2007	
A	27.03.2007		PP	27.03.2007	JF
Rev	Crea. date dd-mm-yy	Revision description	Draw. Init	Appr. date dd-mm-yy	Appr. Init
F	<p>Struers Pederstrupvej 84 DK-2750 Ballerup Copenhagen Denmark Phone : +45 44600 800 Fax : +45 44600 804</p>	Material:	Scale: 1:1	Format: A4	Tolerance: DS/ISO 2768 - Surface treat.:
		ID: 15920014	Description: Valve - assembly	Rev: B	

B	23-06-08	Reference letters added	BRY		
A	21-05-07		PP	21-05-07	JF
Rev	Crea. date dd-mm-yy	Revision description	Draw. Init	Appr. date dd-mm-yy	Appr. Init
			Scale: 1:1	Format: A4	Tolerance: DS/ISO 2768 - Surface treat.:
ID:		Description: 15922000 Air diagram Citovac - ejector			Rev: B

Struers

Pedestrupvej 84
DK-2750 Ballerup
Denmark
Phone: +45 44600 800
Fax: +45 44600 804

1

2

3

4

A


B

C




D

E


F

B	23-06-08	Reference letters added	BRY		
A	21-05-07		PP	21-05-07	JF
Rev	Crea. date dd-mm-yy	Revision description	Draw. Init	Appr. date dd-mm-yy	Appr. Init
		Material:	Scale: 1:1	Format: A4	Tolerance: DS/ISO 2768 - Surface treat.:
Pederstrupvej 84 DK-2750 Ballerup Copenhagen Denmark Phone: +45 44600 800 Fax: +45 44600 804		ID: Description: 15922001 Air diagram Citovac - no ejector			Rev: B







Note #1:
X20 is only present in the variant with external vacuum pump

REV A: OH (2007/07/02)	STRUERS A/S Pederstrupvej 84 DK-2750 Ballerup Denmark		
	CitoVac, BLOCK DIAGRAM (FOR VARIANT WITH EJECTOR AND EXTERNAL VACUUM PUMP)		
Size A3	CAGE Code	DWG NO 15923050	Rev A
Monday, July 02, 2007	Scale	OH	Sheet 1 of 1

Main Supply Input:
100-120V, 50/60Hz
200-240V, 50/60Hz

COLOR CODE:
BK = black
BN = brown
RD = red
OG = orange
YE = yellow
GN = green
BU = blue
VT = violet
GY = grey
WH = white
PK = pink
BE = beige
RO = rose
GNYE = yellow-green







Note #1:
W2 and X20 is only present in the variant
with external vacuum pump

Note #2:
X4-2 and X4-4 connected - variant with external vacuum pump
X4-2 and X4-4 not connected: variant with the ejector



REV A: OH (2007/07/02)	STRUERS A/S Pædelstrøvej 64 DK-2750 Ballerup Denmark		
CitoVac, CIRCUIT DIAGRAM (FOR VARIANT WITH EJECTOR AND EXTERNAL VACUUM PUMP)			
Size A2	CAGE Code	DWG NO 15923100	Rev I A
Scale	OH	Sheet 1	of 1

Friday, August 24, 2007

Bottom Plate (view from inside of the machine)

Front panel (view from the bottom side)

- COLOR CODE:
- BK = black
 - BN = brown
 - RD = red
 - OG = orange
 - YE = yellow
 - GN = green
 - BU = blue
 - VT = violet
 - GY = grey
 - WH = white
 - PK = pink
 - BE = beige
 - RO = rose
 - GNYE = yellow-green

Note #1:
W2 and X20 is only present in the variant with external vacuum pump

REV A: OH (2007/07/02) REV B: AKR (23.06.08) Names on A1 (L,N) changed		STRUERS A/S Føderstrupvej 84 DK-2750 Ballerup Denmark	
CitoVac, WIRING DIAGRAM (FOR VARIANT WITH EJECTOR AND EXTERNAL VACUUM PUMP)			
Size A2	CAGE Code	DWG NO 15923450	Rev 1 B
Monday, June 23, 2008	Scale	OH	Sheet 1 of 1

English

Declaration of Conformity**Manufacturer,
responsible for
Technical File**Struers A/S
Pederstrupvej 84
DK-2750 Ballerup, Denmark
Telephone +45 44 600 800

Herewith declares that

<i>Product Name:</i>	CitoVac
<i>Type No:</i>	592
<i>Machine Type:</i>	Vacuum impregnation apparatus

is in conformity with the provisions of the following directives:

Safety of Machinery 2006/42/EC according to the following standard(s):
EN ISO 12100:2011, EN ISO 13849-2:2014, EN 60204-1:2006/AC:2010, EN 61010:2010,
EN 13218:2002+A1:2008/AC:2010.**EMC-Directive** 2004/108/EC according to the following standard(s):
EN 61000-6-2:2005, EN 61326-1:2013.**RoHS** 2011/65/EU according to the following standard(s):
EN 50581:2012.**Supplementary
Information** The equipment complies with the following standards:
UL508, NFPA70:2014; NFPA79:2012, FCC part 15, subpart B.**The above has been declared according to the global method, module A**

Date: 13.02.2015

Christian Skjold Heyde,
Vice President, R & D and Production, Struers A/S

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Stellvertretender Geschäftsführer, Entwicklung und Produktion, Struers A/S
Technique** Struers A/S
Vice- President, R & D et Production, Struers A/S